

Magnetic tapes supported by Advanced Logic Solutions

The following chart lists the magnetic tape media and formats that Advanced Logic Solutions can process. Tapes should be initialized to one of the densities shown in the chart. Under most situations, the lowest capacity tape that your tape drive can write to will be sufficient to hold the requested data files library.

1. You must provide a “New volume identifier” (NEWVOL) on the “Initialize tape” (INZTAP) command.
2. When saving a library to the initialized tape, the “Target release” (TGTRLS) parameter on the “Save Library” (SAVLIB) command must be V5R4M0 or earlier.

Format (Density)	Capacity	Media	Cartridge
			1/4” SLR60
*MLR3	25 GB	MLR3-25GB	RW
*QIC2DC	5 GB	DC9250	R
*QIC2GB	2.5 GB	DC9250	R
*QIC4DC	8 GB	SLR5-4GB	R
*QIC4GB	4 GB	SLR5-4GB	R
*QIC5010	2 GBB	MLR1-2GB	R/W
*QIC5010	16 GB	MLR1-16GB	R/W
*QIC5010	13 GB	DC5010	R/W
*SLR100	5 GB	SLR100-5GB	RW
*SLR60	30 GB	SLR60-30GB	RW